


Instructions for Use

III BODY COMPOSITION ANALYZER

MC-780MA-N



Separate Version

Column Mounted Version

<Usage Conditions>

| | |
|----------------------------|------------------------------------|
| Temperature | : 5-35°C |
| Relative Humidity | : 30-80% (without condensation) |
| Maximum Altitude | : 2,000m ASL |
| Atmospheric Pressure Range | : 86-106kPa |

<Storage Conditions>

| | |
|----------------------------|---------------------------------|
| Temperature | : -10-50°C |
| Relative Humidity | : 10-90% (without condensation) |
| Atmospheric Pressure Range | : 70-106kPa |

To avoid malfunctions, avoid storing the equipment in a place with direct sunlight, significant temperature changes, a risk of dampness, a large amount of dust or a risk of vibration or impact, or in the vicinity of flames.



Please read this Instruction Manual carefully and keep it for future reference.

Contents

en

Before Use

| | |
|--|----------|
| Part Names & Connection Procedure | 3 |
| Power Supply | 5 |
| Settings | 6 |

How to Use

| | |
|--|-----------|
| Taking a Measurement | 8 |
| Body Composition Mode | 8 |
| Measurement Results | 10 |
| Scale Mode | 12 |
| Weight Lock Mode | 13 |
| General Instructions for Body Composition Measurement | 14 |
| Output and Storage of Measurement Results | 15 |


If Necessary

| | |
|------------------------------|-----------|
| Troubleshooting | 17 |
| Specifications | 18 |

Part Names & Connection Procedure


Symbols and their Meanings

| | | | |
|--|--|--|---|
| | Conformity with Medical Device Directive 93/42/EEC | | Serial interface |
| | Alternating current | | Direct current |
| | Input, Output | | Class II Equipment |
| | Manufacturer | | Type BF applied part - Grips and platform |
| | Caution Refer to the attached notes. | | WEEE - Waste Electrical and Electronic Equipment Directives |
| | SD card | | For indoor use only |
| | Serial number | | See the instructions |




Checking the level

- * For accurate measurement, place the machine as level as possible.
- * Rotate the adjustable feet in 4 positions for adjustment so that the bubbles of the level gauge reach the centre.




Status when the level gauge is viewed from above


Control unit Connection of Plugs



Platform Connection of Plugs



Display and Keys



Meanings of the LED Indicators and Keys

| | | | |
|--|-----------------------------------|---|---|
| | Turn ON / OFF the power | | Display Body Fat (percentage and mass) *Not measured value but calculated value |
| | Printout result | | Display Body Muscle (percentage and mass) *Not measured value but calculated value |
| | Set various functions | | Display Body Water (percentage and mass) *Not measured value but calculated value |
| | Reset zero point | | Display Visceral Fat Rating |
| | Set preset value (Clothes weight) | | Display Basal Metabolic Rate |
| | Select measurement mode | <input type="radio"/> | Indicates whether "Standard mode or Athletic mode" is selected as the body type. |
| | | <input type="radio"/> <input type="radio"/> | Indicates whether "Male or Female" is selected as the gender |
| | | <input type="radio"/> 5 - 99 | Enter the age between "5-99 years old" |
| | | <input type="radio"/> 90.0 - 249.9cm | Enter the height between "90.0-249.9 cm" |
| | | | Confirms the entered numerical value. |


Power Supply

Turning the Main Power ON/OFF

Turning the main power ON.


Press the  key to turn on the power.

The initial screen is displayed.




Turning the main power OFF.

Press the  key to turn the power off.



Emergency Shut Down

Keep the area around the plug socket clear during operation of the equipment in case of an emergency.




Settings

en

Before Use

1 Press the  key to change the mode.

The setting screen is displayed.



Note

• The  key cannot be used when the scale is measuring weight or impedance.

2 Select the setting item from the list below.
Enter numerical values and Press the  key.

Setting item List  → Save changes and return  → Return without saving


| No. | Setting Item | Default | Page | No. | Setting Item | Default | Page |
|-----|---|---------|------|-----|---|---------|------|
| 0 | Check the software version | - | - | 10 | Height input unit 0: 0.1cm increments 1: 1cm increments | 0 | P.9 |
| 1 | Date and time | - | P.6 | 11 | Automatic determination time of input information 0-9 second(s) * 0: disables this function | 0 | - |
| 2 | Number to be printed automatically (body composition mode) 0-2 sheet(s) | 0 | P.16 | 18 | Target body fat ratio input 0: off 1: on | 0 | P.9 |
| 3 | Number to be printed automatically (scale mode and weight lock mode) 0-2 sheet(s) | 0 | P.16 | 19 | Printout language 1: English, 2: French, 3: German, 4: Italian, 5: Spanish | 1 | - |
| 4 | Select printing paper 0: pre printed paper 1: plain paper | 1 | - | 20 | Adjust printing position 2: ↓, 4: ←, 6: →, 8: ↑ | - | - |
| 5 | Beep sound 0: off 1: on | 1 | - | 21 | Timeout function of result display 0: disable 1: enable | 0 | - |
| 6 | ID number (manual input) 0: off 1: on | 0 | P.8 | 30 | Data output format to PC 0: default format 1: BC-418 format | 0 | P.15 |
| 7 | ID number (automatic count up) 0: off 1: on | 0 | P.8 | 31 | Display direction 0: default 1: reversal position | 0 | - |
| 8 | Measurement flow 0: two step flow (measure body weight first) 1: one step flow (enter personal information first) | 0 | P.8 | 32 | BMI desirable range 0: 18.5-24.9 1: 18.5-22.9 2: 18.5-23.9 | 0 | - |
| 9 | Body type selection (athletic mode) 0: off 1: on | 0 | P.8 | 60 | Reading stored measurement result | - | P.16 |

Note

Enter the year, month, day, hour and minute.
The date format is "yyyy mm dd hh:mm"

Example 6:47 pm, 24th September, 2024
"2024" "09 24" "18:47"


To enter a number with 1 digit (0-9), press "0" first.



Select the Measurement Mode


Select the measurement mode by pressing the  key.

The measurement mode is switched in the following order when the  key is pressed:




The "" mark is displayed when the scale mode is selected.


Body Composition Mode



Scale Mode



Scale Mode (Weight Lock Mode)



Taking a Measurement

en

How to Use



Body Composition Mode (Two step flow Measure body weight first)

1 Enter the preset tare value (clothes weight) and press **Enter** key

Turn the power on, and check that the "PT" mark flashes.

Enter the preset tare value (clothes weight)



The tare value range is 0.0–10.0kg.



2 Measure body weight

Step onto the scale, after "STEP on" is displayed.

"NET" is displayed when you have entered a tare value (clothes weight).



Note

When the One step flow is selected, enter personal information first. → P.6 Setting 8


If the scale does not detect a load, press **Clothes Weight (PT)** to switch to "Input tare value".

If the scale detects the load, press **Clothes Weight (PT)** to display the entered tare value.

3 Enter ID number

If the ID number function is set to OFF, the scale switches directly to "Select body type" without requesting an ID number.


The ID number range : 0–9999999999999999



4 Select body type **Standard** **Athletic**

If the body type selection (Athletic mode) is set to OFF, the scale switches directly to "Select gender" without requesting a body type.

The "⚡" mark is displayed when the Athletic mode is selected.




Note

Athletic mode setting → P.6 Setting 9



Athlete condition → P.14

Body Composition Mode

5 Select gender





6 Enter age and press Enter key



The age range : 5-99

7 Enter height and press Enter key



The height range : 90.0-249.9 cm


8 Set target body fat ratio

If the target body fat ratio function is set to OFF, the scale switches directly to "Measuring body composition" without entering the target value.

The body fat ratio range: 4-55 %

Note

Target Body Fat Ratio setting → P.6 Setting 18






9 Measuring body composition

Hold the grip, when "grip on" is displayed.

The scale starts measuring impedance after you have entered all of the personal information.

The segmental impedance values are measured in order.

The scale displays the measurement results after measuring the whole body impedance and segmental impedance.




Taking a Measurement (continued)

en

How to Use

Measurement Results (Body Composition Mode)

Output Measurement Results (Body Composition Mode)



The measurement results are displayed on the LCD after measurement is completed. The results are output to the PC immediately after measurement is completed.

Press the key to select the measurement display.



The measurement display is switched in the following order by pressing the key.

"Body Fat" → "Muscle" → "Body Water" → "Visceral fat rating" → "Basal metabolic rate" → "BMI"

Body Fat

Fat icon

Segmental fat %



Body weight (kg)

Body fat %

Body fat mass (kg)



Evaluation of your results compared to average values.

Fat rating

Muscle

Muscle icon

Segmental Muscle Mass



Body weight (kg)

Body muscle %


Body muscle mass (kg)

Evaluation of your results compared to average values.


Muscle mass rating

Measurement Results (Body Composition Mode)


Body Water




Visceral Fat Rating



Basal Metabolic Rate (BMR)



Body Mass Index (BMI)



Note

“○” flashes if measuring posture is incorrect.
The results can only be used as reference data.



Taking a Measurement (continued)

en

Scale Mode

How to Use

The "I" icon is displayed when the scale mode is selected.

1 Enter the ID number


If the ID number function is set to OFF, the scale starts measuring weight immediately.

Enter the ID number.

The ID number range is from 0-9999999999999999.

Press **PT** to enter tare value (clothes weight).



The tare value range is 0.0-10.0kg.



2 Measure body weight


"NET" is displayed when you have entered a tare value (clothes weight).

Step onto the scale, after "STEP ON" flashes.



The "Stabilised" icon (●) appears when the load is stable. The "Stabilised" icon (○) disappears when the load is unstable.

When the load is stable, the measurement results are output via the USB port and stored on an SD card.




Weight Lock Mode

The "🔒" icon and "Weight Lock Mode" are displayed when the weight lock mode is selected.

Enter the preset tare value (clothes weight)


Enter the preset tare value (clothes weight)
The tare value range is 0.0-10.0kg.



Enter the ID number

If the ID number function is set to OFF, the scale starts measuring weight immediately.



Enter the ID number.
The ID number range is from 0-9999999999999999.



Measure body weight

"NET" is displayed when you have entered a tare value (clothes weight).

Step onto the scale, after
"StEP on" flashes .



The "Stabilised" icon (●) appears when the load is stable.
The "Stabilised" icon (○) disappears when the load is unstable.

When the load is stable, the measurement results are output via the USB port and stored on an SD card.

Taking a Measurement (continued)

en

How to Use

General Instructions for Body Composition Measurement

Athletic Mode

- **Recommended for those who are 18 years or older and meet the following conditions.**
 - People who carry out 12 hours or more of cardiovascular exercise a week.
 - People who belong to a sport team or a sport organization with the aim of participation in competition, etc.
 - People who are professional athletes.

Target Body Fat

- **A target body fat % should be set by a professional only. Tanita is not responsible for setting the appropriate target body fat % for specific individuals.**

Attention

- **Posture when measuring**
 - Stand with both feet parallel on the electrodes.
 - Stand without bending knees.
- **The age input range is 5–99 years old.**
Input age 99 for those who are 100 years or older.
- **Contact condition**
 - If your inner thighs are touching or your arms make skin-to-skin contact, the results can only be used as reference data.

Note

- Inaccurate results may be reported after excessive food/fluid intake, or after periods of intense exercise.
- If clothes weight is input, clothes weight is subtracted from the weight measurements.

Output and Storage of Measurement Results



Data Output

Outputting measurement results via the USB or RS-232C port

The results are output to the PC immediately after measurement is completed.

Data is output in CSV format.

- USB-B port is located on the back of the control box.
- *communication cable is not included in the product.
- You must install the necessary driver onto your PC, available download from <https://www.tanita.eu/>



Data Storage

Data storage of measurement results

The results are saved to the SD card immediately after measurement is completed.

The  icon is displayed when a valid SD card is installed.


A new file is created on the SD card for each day.

The file is created using the measurement date and time as the file name, as shown below. "YYYYMMDD" (year, month and date)

If "error 12" is shown in the display, this indicates that there is not enough free space left on the SD card.

You can store approximately 10,000 measurement results (weight only mode) on a 2GB SD card.

- Do not remove the SD card when it is storing or reading data.
- Do not turn off the main power when the SD card is storing or reading data.



Printing Data

Printing Measurement Results

The results are printed from the printer immediately after measurement is completed.

The number of sheets that are printed is set in  .

Press the  key to print additional copies.

You can print out the result if  is lighted.

Output and Storage of Measurement Results (continued)

en


How to Use

Reading Stored Measurement Result Data

Reading stored measurement result data (→P.6 Setting 60)

Select the file that you want to load using the  or  key and  or  key and press **Enter**.

Displaying stored data
Date and time when the displayed data was measured
Year, Month and date
Hour, Minute and Second




Switch the display by pressing the  key.

 → Print

Enter → Next data

CE → Exit



Troubleshooting

Please check the following before requesting repair.

| | Problem | Solution |
|-------------|---|---|
| Measurement | Error with impedance measurement | <ul style="list-style-type: none"> • Hold the handgrips and electrodes firmly with fingers and palms. • Stand barefoot on the platform for measurement. • If the soles of the feet are dry, use a dropper to apply about 0.5 mL of water before measurement. • Check the input details. |
| | Error with zeroing | <ul style="list-style-type: none"> • Turn off the power and remove anything on the platform, then turn on the power and try measuring again. |
| | The weight value does not stabilize. | <ul style="list-style-type: none"> • Is the equipment placed on a vibrating surface? • Is the measuring platform tilted? • Is something blocking the gap in the measuring platform? • Remove any inserted objects. |
| Display | Nothing is displayed, even when the power is turned on. | <ul style="list-style-type: none"> • Confirm that the power is connected correctly. |
| | "-----" is displayed. | <ul style="list-style-type: none"> • The measured weight exceeds weighing capacity. |
| | "○" flashes. | <ul style="list-style-type: none"> • "○" flashes if measuring posture is incorrect. The results can only be used as reference data. |

If Necessary

Specifications

en

If Necessary

| | | |
|------------------------------|--|---|
| Model Number | | MC-780MA-N |
| Classification | | MDD: CLASS IIa NAWI: CLASS III |
| Power Source | | 100–240V AC |
| Rated Power | | 36VA |
| Impedance Measurement | Measurement System | Multi-frequency 8 electrodes |
| | Measurement Frequency | 5kHz/50kHz/250kHz |
| | Measurement Current | 90µA or less |
| | Electrode Materials | Feet: Stainless steel/Handgrips: plated |
| | Measurement Part | Whole body/Right arm/Left arm/Right leg/Left leg |
| | Measurement Range | 75.0–1,500.0Ω (0.1Ω increments) |
| | Accuracy at First Calibration | ± 2% |
| Weight Measurement | Measurement System | Strain gauge load cell |
| | Maximum Capacity | 270kg (including preset tare value) |
| | Minimum Graduation | 0.1kg |
| | Accuracy at First Calibration | ± 0.2kg |
| Display | | Dual LCD screen |
| Interface | | USB Type-B port (device) RS-232C USB Type-A port (for USB printer) SD card |
| Product Weight | | Column mounted version: 15.5kg/Separate version: 11.1kg |
| Product Size | Platform | Platform size: 360 x 360mm Height: 94mm |
| | Height (Column Mounted Version) | 1165mm |
| Time accuracy | | Within ±1 minute per month (+18–+28°C) |

| | | | | |
|----------------------------|---------------------------------------|---|---------------------------------|-----------------|
| Input items | Clothes Weight | | 0–10.0kg (0.1kg increments) | |
| | Serial No. | | Within 16 digits | |
| | Gender | | Female/Male | |
| | Body Type | | Standard/Athletic* ¹ | |
| | Age | | 5–99 years old | |
| | Height | | 90.0–249.9cm (0.1cm increments) | |
| | Target Body fat % | | 4–55% (1% increments) | |
| Output items* ³ | ID | | 16 digits | |
| | Gender | | Female/Male | |
| | Body Type | | Standard/Athletic* ¹ | |
| | Age | | 5–99 years old | |
| | Height | | 90.0–249.9cm (0.1cm increments) | |
| | Clothes Weight | | 0–10.0kg (0.1kg increments) | |
| | Whole Body Analysis | Weight | 2.0–270.0kg (0.1kg increments) | 5–99 years old |
| | | Fat % | 3.0–75.0% (0.1% increments) | |
| | | Fat Mass | 0.1kg increments | |
| | | FFM | 0.1kg increments | |
| | | Fat Rating | -4–+4 (1 increments) | |
| | | Muscle Mass | 0.1kg increments | 18–99 years old |
| | | Muscle Mass Rating | -4–+4 (1 increments) | |
| | | BMI | 0.1 increments | 5–99 years old |
| | | Bone Mass | 0.1kg increments | 18–99 years old |
| | | Metabolic Age | 1 increments | |
| | | Basal Metabolic Rate | 1kcal/1kJ increments | 5–99 years old |
| | | Visceral Fat Rating | 1–59 (1 increments) | 18–99 years old |
| | | TBW | 0.1kg increments | 5–99 years old |
| | | TBW % | 0.1% increments | |
| | | ECW | 0.1kg increments | 18–99 years old |
| | ICW | 0.1kg increments | | |
| | ECW / TBW | 0.1% increments | | |
| | Segmental Analysis | Muscle Mass | 0.1kg increments | 5–99 years old |
| | | Muscle Mass Rating | -4–+4 (1 increments) | 18–99 years old |
| | | Fat % | 0.1% increments | 5–99 years old |
| | | Fat Mass | 0.1kg increments | |
| | | FFM | 0.1kg increments | |
| | Fat Rating | -4–+4 (1 increments) | 18–99 years old | |
| | Body Balance Evaluation | Physique Rating | | |
| Muscle Mass Balance | | | | |
| Leg Muscle Score | | | | |
| Others | MM/H² | Whole body Muscle Mass (kg)/Height (m) ² | 5–99 years old | |
| | MM/BW | Whole body Muscle Mass (kg)/Body Weight (kg) | | |
| | SMI* | Total muscle mass value of Right arm, Left arm, Right leg and Left leg (kg)/Height (m) ² | | |
| | ASM**/BW | Total muscle mass value of Right arm, Left arm, Right leg and Left leg (kg)/Body Weight (kg) | | |
| | Bioelectrical data | Reactance/Resistance/Phase Angle | | |
| | Contact Condition*² | | | |

*¹ Athletic mode can be selected only for ages 18–99.

*² The results can only be used as reference data. → P.11

*SMI: Skeletal Muscle Mass Index

**ASM: Appendicular Skeletal Muscle Mass

Note

Compatible with SD and SDHC memory cards.
Not compatible with SDXC memory cards.



- This product meets the following requirements;
1. Medical Device Directive (93/42/EEC)
 2. Non-Automatic Weighing Instruments (2014/31/EU)
 3. RoHS Directive (2011/65/EU)

Disposal



This product is an electronic device. Please dispose of this product appropriately, not as general household waste. Be sure to follow the regulations in your area when disposing of this product.

- The product design and specifications may be changed at any time without prior notice.
- SD, SDHC and SDXC Logos are trademarks or registered trademarks of SD-3C LLC in the United States, other countries or both.

<Manufacturer>
TANITA Corporation
1-14-2 Maeno-cho, Itabashi-ku,
Tokyo 174-8630 Japan
TEL: +81-(0)3-3968-7048
www.tanita.co.jp

<EU Representative>
EU Importer
TANITA Europe B.V.
Hoogoorddreef 56-E, 1101 BE
Amsterdam, the Netherlands
TEL: +31-(0)20-560-2970
FAX: +31-(0)20-560-2988
www.tanita.eu

TANITA Corporation of America, Inc.
2625 South Clearbrook Drive,
Arlington Heights, Illinois 60005
U.S.A.
TEL: +1-847-640-9241
FAX: +1-847-640-9261
www.tanita.com

TANITA Health Equipment H.K. Ltd.
Unit 301-303, Wing On Plaza, 3/F., 62 Mody
Road,
Tsimshatsui East, Kowloon, Hong Kong
TEL: +852-2834-3917
FAX: +852-2838-8667
www.tanita.asia

TANITA (Shanghai) Trading Co., Ltd.
Room 8005, 877 Huai Hai Zhong Lu,
Shanghai,
The People's Republic of China
TEL: +86-21-6474-6803
FAX: +86-21-6474-7901
www.tanita.com.cn

<UK Representative>
UK Importer
TANITA Europe B.V. Manchester Office
111 Piccadilly, Manchester, M1 2HY, UK
TEL: +44-161-6380926
www.tanita.eu